


the new *force* in laser therapy



lightforce™

LightForce™ EX Therapy Laser

with Influence™ Technology and the Empower™ Delivery System

EXpert Design

The LightForce™ EX therapy laser by LiteCure® Medical is designed to help you achieve more consistent positive outcomes for your most difficult clinical conditions. The LightForce™ EX therapy laser is the perfect combination of power, intelligence and versatility.



get your patients back in *action*



Influence™ Technology

maximizing outcomes

Delivery of the optimal therapeutic dose to the damaged tissue is the key to effective results.

Muscle tone, size and pigmentation significantly affect the penetration of all therapy lasers.

Influence™ technology provides patient specific protocols based on individual characteristics. Personalize your treatments with Influence™ and achieve consistent results.

unit operation

current treatment

h: large
b: medium
s: medium
bursitis
hip

lightforce

Energy Delivered (J) 6840
Total Treatment Energy (J) 8640



Time 12:00
Power (W) 12
Pulsed Rate (Hz) CW

LightForce™ Power

feel the difference

Featuring up to 15W of therapeutic power, the LightForce™ EX therapy laser facilitates fast and efficient treatments deep in targeted tissues. Experience the soothing warmth and powerful benefits of the LightForce™ EX therapy laser.

EMPOWER Delivery System™

Empower DS™ provides the versatility you need to treat your toughest cases. Four available treatment heads, the patented laser-massage ball and the toughest fiber available make the LightForce™ advantage yours.

the advantage is yours

lightforce™

EX
expert

a LiteCure® Medical therapy laser



Specifications

Laser Type	Class IV, Solid State
Laser Wavelength	980/810 nm
Laser Power	0.5W - 15W
Operating Modes	CW or Pulsed
Aiming Beam	650 nm, 4.5 mW
Dimensions	15.5" x 9.75" x 10.5"
Weight	17 lbs.

Customer Support

LiteCure Medical boasts the industry's leading customer support programs:

Clinical Resources. When you have clinical questions we have qualified clinicians and specialists on call to provide you with answers.

Educational Resources. Third party educational resource, the American Institute of Medical Laser Applications provides foundational knowledge.

Marketing Support. LiteCure Medical provides you with the marketing resources you need to successfully implement a laser therapy program in your practice.

Loaner Program. Our exceptional loaner program ensures that you never spend a day without your most versatile treatment modality.

Warranty

LiteCure's superior platinum warranty keeps your investment safe for 24 months. The comprehensive warranty includes coverage for the therapy unit, hand pieces, and fiber.

P | 302.709.0408
F | 302.709.0409



litecuremedical.com | info@litecure.com | 250 Corporate Blvd | Suite B | Newark, DE 19702